
Lithium sulphate monohydrate, Hi-LR™

GRM567


Product Identifier

| | | |
|-------------------|---|---------------------------------------------------|
| CAS No. | : | 10102-25-7 |
| EC No. | : | 233-820-4 |
| Molecular Formula | : | Li ₂ SO ₄ .H ₂ O |
| Molecular Weight | : | 127.96 |
| HS Code | : | 2833 29 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|--------------------------------|---|---------------------------------|
| Appearance | : | White crystals or powder |
| Solubility | : | 100 mg soluble in 1 mL of water |
| Assay (Ion-Exchange Titration) | : | 98.00 - 102.00% |

Safety Information

| | | |
|---------------------|---|------------------------------------------------------------------------------------|
| Hazard Pictogram(s) | : |  |
| Signal Word | : | Warning |
| Hazard Statement(s) | : | H302 |
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | OJ6419000 |
| WGK | : | 1 |