
Chromium carbonate, Hi-LR[™]

GRM5742

Product Identifier

CAS No.	:	29689-14-3
EC No.	:	249-781-1
Molecular Formula	:	Cr ₂ O ₃ .xCO ₂ .yH ₂ O
HS Code	:	2836 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Blue-green or light green crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of hot hydrochloric acid
Assay (Cr, Iodometry/ICP)	:	50.00 - 55.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1