
Ethyl benzoate

GRM5887

Product Identifier

CAS No.	:	93-89-0
EC No.	:	202-284-3
Molecular Formula	:	C ₉ H ₁₀ O ₂
Molecular Weight	:	150.17
HS Code	:	2916 31 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to pale yellow liquid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n _{20/D})	:	1.504 - 1.507
Density (at 25°C)	:	1.040 - 1.050 g/mL
Assay (GC)	:	min. 99.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH0200000
WGK	:	1