
Potassium sodium tartrate tetrahydrate

GRM598

Product Identifier

| | | |
|-------------------|---|-----------------------------------------------------------------------|
| CAS No. | : | 6381-59-5 |
| EC No. | : | 206-156-8 |
| Molecular Formula | : | C ₄ H ₄ KNaO ₆ ·4H ₂ O |
| Molecular Weight | : | 282.22 |
| Synonym | : | Sodium potassium tartrate tetrahydrate; Rochelle salt; Seignette salt |
| HS Code | : | 2918 13 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|----------------------------|---|-------------------------------------------------|
| Appearance | : | Colourless or white crystals or powder or solid |
| Solubility | : | 100 mg soluble in 1 mL of water |
| pH (5% in water at 20°C) | : | 5.50 - 8.50 |
| Chloride (Cl) | : | <= 0.005% |
| Iron (Fe) | : | <= 0.001% |
| Sulfate (SO ₄) | : | <= 0.005% |
| Assay (NT) | : | 98.00 - 102.00% |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | Not available |
| WGK | : | 1 |