

GRM5989**Iodoform****Product Number****Packing**

GRM5989 : 250G

Product Information

Product Code : GRM5989
Product Name : Iodoform
Synonym : Triiodomethane
Molecular Formula : CHI_3
Molecular Weight : 393.73
CAS No. : 75-47-8
EC No. : 200-874-5
HS Code : 2903 39 90
Shelf Life : 4 years

Technical Specification


Appearance : Bright yellow to greenish yellow crystals or powder
Solubility : 33.3 mg soluble in 1 mL of ethanol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 115 - 125°C
Assay (Argentometric-Back Titration/GC) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302+H312+H332-H315-H319-H335
Precautionary Statement(s) : P261-P280-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 20/21/22-36/37/38
S-Phrase(s) : 26-36/37
WGK : 3
RTECS : PB7000000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods



Product Information

IATA

: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein