

GRM6048 Cholesterol, Hi-LR™

Product Number Packing

GRM6048 : 100G
GRM6048 : 500G

Product Information

Product Code : GRM6048
Product Name : Cholesterol, Hi-LR™
Synonym : 5- Cholesten-3β-ol
Molecular Formula : C₂₇H₄₆O
Molecular Weight : 386.65
CAS No. : 57-88-5
EC No. : 200-353-2
HS Code : 2906 13 10
Shelf Life : 3 years

Technical Specification

Appearance : White to off-white crystals or powder or flakes
Solubility : 100 mg soluble in 1 mL of chloroform
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 142 - 150°C
Assay (GC) : 95.00 - 100.00%

Risk and Safety Information

WGK : 1
RTECS : FZ8400000
Storage Temperature(°C) : Store at 2 - 8°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods