

Methyl benzoate

GRM6053

Product Identifier

CAS No.	:	93-58-3
EC No.	:	202-259-7
Molecular Formula	:	C ₈ H ₈ O ₂
Molecular Weight	:	136.15
HS Code	:	2916 31 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to pale yellow liquid
Solubility	:	1 mL miscible in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n _{20/D})	:	1.515 - 1.518
Density (at 25°C)	:	1.083 - 1.093 g/mL
Assay (GC/NT)	:	99.00 - 100.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H302
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DH3850000
WGK	:	1