

GRM6096 α -Naphthylamine, Hi-AR™**Product Number Packing**

GRM6096 : 100G

Product Information

Product Code : GRM6096
Product Name : α -Naphthylamine, Hi-AR™
Synonym : 1-Naphthylamine; 1-Aminonaphthalene
Molecular Formula : C₁₀H₉N
Molecular Weight : 143.19
CAS No. : 134-32-7
EC No. : 205-138-7
HS Code : 2921 45 00
Shelf Life : 4 years

Technical Specification

Appearance : Tan to purple crystals or powder or solid
Solubility : 33.3 mg soluble in 1 mL of methanol
FTIR (KBr disc) : Matches with the standard pattern
Melting range : 47 - 51°C
Assay (GC/NT) : 99.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H302-H350-H411
Precautionary Statement(s) : P201-P273-P308+P313
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 45-22-51/53
S-Phrase(s) : 53-45-61
WGK : 3
RTECS : QM1400000
Flash Point(°F) : 314.6 °F
Flash Point(°C) : 157 °C
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**



Product Information

UN No.	:	2077
Class	:	6.1
Packaging Group	:	3
Marine Pollutant	:	Yes
ADR/RID	:	2077 6.1/PG 3
IMDG	:	2077 6.1/PG 3
IATA	:	2077 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein