

Sodium desoxycholate

GRM6216


Product Identifier

CAS No.	:	302-95-4
EC No.	:	206-132-7
Molecular Formula	:	C ₂₄ H ₃₉ NaO ₄
Molecular Weight	:	414.55
HS Code	:	2918 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to cream crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (2% in water at 20°C)	:	7.50 - 9.50
FTIR	:	Matches with the standard pattern
Specific rotation [alpha 20/D]	:	+42.0° to +46.0° (c = 2% in water, at 20°C)
Loss on drying (at 105°C, 4 hr)	:	<= 5.00%
Assay (NT, on dry basis)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302
Precautionary Statement(s)	:	P264- P270- P301 + P312- P330- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	FZ2250000
WGK	:	1