

RM6224 Sodium persulphate**Product Number Packing**

RM6224 : 500G

Product Information

Product Code : RM6224
Product Name : Sodium persulphate
Synonym : Sodium peroxodisulfate
Molecular Formula : $\text{Na}_2\text{S}_2\text{O}_8$
Molecular Weight : 238.10
CAS No. : 7775-27-1
EC No. : 231-892-1
HS Code : 2833 40 00
Shelf Life : 4 years

Technical Specification

Appearance : White crystals or powder
Solubility : 100 mg soluble in 1 mL of water
Assay (KMnO₄ Titration/Iodometry) : 99.00 - 102.00 %

GHS Safety Information

Hazard Statement(s) : H272-H302-H315-H317-H319-H334-H335
Precautionary Statement(s) : P220-P261-P280-P305+P351+P338-P342+P311
Signal Word : Danger
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 8-22-42/43
S-Phrase(s) : 22-36/37-45
WGK : 1
RTECS : SE0525000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

UN No. : 1505
Class : 5.1
Packaging Group : 3
Marine Pollutant : No
ADR/RID : 1505 5.1/PG 3



Product Information

IMDG : 1505 5.1/PG 3
IATA : 1505 5.1/PG 3

A-406, Bhaveshwar Plaza, LBS Marg, Mumbai-400 086, India. Phone:+91-022-4095 1919 Fax:+91-022-4095 1920 Gram: STERILITY.
email:info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein