

**GRM6226 Sodium phytate, From Rice**

**Product Number Packing**

---

GRM6226 : 25G  
GRM6226 : 100G

**Product Information**

---

Product Code : GRM6226  
Product Name : Sodium phytate, From Rice  
Synonym : Phytic acid sodium salt  
Molecular Formula :  $C_6H_6Na_{12}O_{24}P_6$   
Molecular Weight : 923.8  
CAS No. : 14306-25-3  
EC No. : 238-242-6  
Shelf Life : 4 years

**Technical Specification**

---

Appearance : White to off-white crystals or powder  
Solubility : 33.3 mg soluble in 1 mL of water

**Risk and Safety Information**

---

WGK : 3  
Storage Temperature(°C) : Store below 30°C

**Transport Information**

---

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods