

---

**di-Sodium tartrate dihydrate, Hi-AR™/ACS****GRM6231****Product Identifier**

---

CAS No.	:	6106-24-7
EC No.	:	212-773-3
Molecular Formula	:	C <sub>4</sub> H <sub>4</sub> Na <sub>2</sub> O <sub>6</sub> .2H <sub>2</sub> O
Molecular Weight	:	230.08
HS Code	:	2918 13 90
Storage	:	Below 30°C
Shelf life	:	4 years

**Technical Specification**

---

Appearance	:	White crystals or powder or needles
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	7.00 - 9.00
Chloride (Cl)	:	<= 0.005%
Iron (Fe)	:	<= 0.001%
Loss on drying (at 150°C, 4 hr)	:	15.61 - 15.71 %
Assay (NT, on dry basis)	:	99.00 - 101.00%

**Safety Information**

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3