

---

## Strontium sulphate, Hi-LR<sup>TM</sup>

**GRM6240**

### Product Identifier

---

CAS No.	:	7759-02-6
EC No.	:	231-850-2
Molecular Formula	:	SrSO <sub>4</sub>
Molecular Weight	:	183.68
HS Code	:	2833 29 50
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	White to light yellow crystals or powder
Solubility	:	10 mg soluble in 10 mL of dilute nitric acid
Assay (EDTA Titration)	:	min. 98.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1