

Succinic anhydride

GRM6243


Product Identifier

CAS No.	:	108-30-5
EC No.	:	203-570-0
Molecular Formula	:	C ₄ H ₄ O ₃
Molecular Weight	:	100.07
HS Code	:	2917 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to almost white crystals or powder flakes or pellet or chips
Solubility	:	33.3 mg soluble in 1 mL of hot methanol
FTIR	:	Matches with the standard pattern
Melting range	:	116 - 121°C
Assay (NaOH Titration/GC)	:	min. 99.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H314- H317- H334
Precautionary Statement(s)	:	P260- P280- P303+P361+P353- P304 + P340- P305+P351+P338- P310
UN No.	:	3261
Class	:	8
Packing Group	:	II
RTECS	:	WN0875000
WGK	:	1