

o-Toluidine, Hi-LR®

GRM6270


Product Identifier

CAS No.	:	95-53-4
EC No.	:	202-429-0
EC Index No	:	612-091-00-X
Molecular Formula	:	C ₇ H ₉ N
Molecular Weight	:	107.15
Synonym	:	2-Methylaniline
HS Code	:	2921 43 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to yellow to reddish brown liquid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n _{20/D})	:	1.571 - 1.573
Density (at 25°C)	:	0.995 - 1.005 g/mL
Assay (GC/NT)	:	99.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H301+H331- H319- H350- H400
Precautionary Statement(s)	:	P261- P280- P301 + P310- P304 + P340- P403 + P233- P501
UN No.	:	1708
Class	:	6.1
Packing Group	:	II
RTECS	:	XU2975000
WGK	:	3