

**GRM6308 Tungstic acid, Hi-AR™**

**Product Number Packing**

---

GRM6308	: 100G
GRM6308	: 500G

**Product Information**

---

Product Code	: GRM6308
Product Name	: Tungstic acid, Hi-AR™
Molecular Formula	: H <sub>2</sub> WO <sub>4</sub>
Molecular Weight	: 249.85
CAS No.	: 7783-03-1
EC No.	: 231-975-2
HS Code	: 2825 90 40
Shelf Life	: 4 years

**Technical Specification**

---

Appearance	: Yellow to greenish-yellow to yellow-orange powder
Solubility	: Insoluble in water
Assay (Titration)	: 99.00 - 102.00%

**Risk and Safety Information**

---

WGK	: 3
Storage Temperature(°C)	: Store below 30°C

**Transport Information**

---

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods