
Oil of castor

GRM6358

Product Identifier

CAS No.	:	8001-79-4
EC No.	:	232-293-8
HS Code	:	1515 30 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to pale yellow viscous liquid
Solubility	:	33.3 mg soluble in 1 mL of alcohol
Refractive index (n _{20/D})	:	1.475 - 1.481
Density (at 25°C)	:	0.952 - 0.966 g/mL
Iodine value	:	83 - 88
Hydroxyl value	:	160 - 168
Saponification value	:	176 - 185

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	NWG