

Oil of eucalyptus

GRM6359

Product Identifier

CAS No.	:	8000-48-4
EC No.	:	283-406-2
HS Code	:	3301 29 24
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to pale yellow liquid
Solubility	:	1 mL miscible in 1 mL of ethanol
Specific rotation	:	-5.00° to +10.00° (c = neat)
Refractive index (n _{20/D})	:	1.4580 - 1.4700
Density (at 25°C)	:	0.897 - 0.924 g/mL
Assay (% Cineole)	:	min. 60.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H226- H304- H315- H317- H411
Precautionary Statement(s)	:	P210- P280- P301 + P310- P331- P370 + P378- P403 + P235
UN No.	:	1993
Class	:	3
Packing Group	:	III
RTECS	:	Not available
WGK	:	2