
Oil of olive

GRM6360

Product Identifier

CAS No.	:	8001-25-0
EC No.	:	232-277-0
HS Code	:	1509 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to green oil
Solubility	:	1 mL miscible in 1 mL of chloroform
Refractive index (n 20/D)	:	1.464 - 1.474
Density (at 25°C)	:	0.905 - 0.915 g/mL
Iodine value	:	79 - 88

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RK4300000
WGK	:	3