

---

## Oil of peppermint

**GRM6361**

---

### Product Identifier

CAS No.	:	8006-90-4
HS Code	:	3301 24 00
Storage	:	Below 30°C
Shelf life	:	4 years


---

### Technical Specification

Appearance	:	Colourless to light yellow or pale greenish-yellow liquid
Solubility	:	1 mL miscible in 3 mL of 70% ethanol
Specific rotation	:	-30.0 to -10.0° (c = neat)
Refractive index (n <sub>20/D</sub> )	:	1.4550 - 1.4670
Density (at 25°C)	:	0.890 - 0.916 g/mL
Assay (% Menthol, GC)	:	30.00 - 55.00%

---

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H317- H411
Precautionary Statement(s)	:	P261- P273- P280- P333 + P313- P362 + P364- P501
UN No.	:	3082
Class	:	9
Packing Group	:	III
RTECS	:	SC6125000
WGK	:	2