
Paraffin liquid, heavy

GRM6362

Product Identifier

CAS No.	:	8012-95-1
EC No.	:	232-384-2
HS Code	:	2710 19 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless, transparent oily liquid
Solubility	:	33.3 mg soluble in 1 mL of ether
Density (at 25°C)	:	0.860 - 0.890 g/mL
Viscosity (at 20°C)	:	110 - 230 mPas (Spindle No. 18 at 12 RPM)

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H304
Precautionary Statement(s)	:	P301 + 310- P331- P405- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	PY8030000
WGK	:	1