
Lecithin soya, 30%

GRM637

Product Identifier

CAS No.	:	8002-43-5
EC No.	:	232-307-2
Synonym	:	L- α -Lecithin, from Soyabean; L- α -Phosphatidylcholine, 30%; Soya lecithin, 30%
HS Code	:	2923 20 10
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to brown powder
Solubility	:	100 mg soluble in 1 mL of chloroform gives clear to slightly hazy solution
Acid value	:	≤ 36.00
Peroxide value	:	≤ 1.00
Heavy metals (as Pb)	:	$\leq 0.002\%$
Water (K.F.)	:	$\leq 1.00\%$
Phosphatides (as acetone insoluble %)	:	min. 95.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OG7565000
WGK	:	1