
Soya lecithin liquid

GRM6398

Product Identifier

CAS No.	:	8002-43-5
EC No.	:	232-307-2
HS Code	:	2923 20 10
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Light brown to brown liquid, having typically soy smell and taste
Water (K.F.)	:	<= 1.0%
Acid value	:	<= 30 mg KOH/gm
Phosphatides (as acetone insoluble %)	:	min. 62.0%
Hexane insoluble	:	<= 0.3%
Peroxide value	:	<= 5 meq/Kg

Safety Information

UN No.	:	Not available
Class	:	-
Packing Group	:	-
RTECS	:	OG7565000
WGK	:	1