
Albumin egg powder

GRM6421

Product Identifier

CAS No.	:	9006-59-1
EC No.	:	232-692-7
HS Code	:	3502 19 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow powder
Solubility	:	33.3 mg soluble in 1 mL of water to give clear to slightly hazy solution
pH (1% in water at 25°C)	:	6.00 - 8.00
Water (K.F.)	:	<= 8.00%
Sulphated ash	:	<= 6.00%
Fat content	:	<= 0.60%
Protein	:	min. 78.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H334
Precautionary Statement(s)	:	P261- P284- P304 + P340- P342 + P311- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AY9384000
WGK	:	1