

**GRM6515 Calcium alginate, from brown algae****Product Number Packing**

GRM6515 : 500G

**Product Information**

Product Code : GRM6515  
Product Name : Calcium alginate, from brown algae  
CAS No. : 9005-35-0  
HS Code : 3913 10 00  
Shelf Life : 4 years

**Technical Specification**

Appearance : White to yellowish-brown to brown granules or powder  
Identification test : 0.5% alkaline sample solution + 2.5% calcium chloride solution gives voluminous, gelatinous precipitate.  
Ash content : 12.00 - 18.00%  
Loss on drying (at 105°C, 4 hr) : <= 20.00%

**Risk and Safety Information**

WGK : 3  
Storage Temperature(°C) : Store below 30°C

**Transport Information**

Marine Pollutant : No  
ADR/RID : Not Dangerous Goods  
IMDG : Not Dangerous Goods  
IATA : Not Dangerous Goods