
Neem seed oil

GRM6541

Product Identifier

CAS No.	:	8002-65-1
HS Code	:	1521 10 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to dark brown or greenish brown viscous oil
Solubility	:	33.3 mg soluble in 1 mL of chloroform
Refractive index (n _{20/D})	:	1.460 - 1.472
Saponification value	:	175 - 210
Iodine value	:	60 - 95

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QO2808000
WGK	:	2