

## D.P.X. mountant (Liquid)

**GRM655**

### Product Identifier

HS Code : 3822 90 90  
Storage : Below 30°C  
Shelf life : 2 years

### Technical Specification

Appearance : Colourless viscous liquid  
Solubility : 1 mL miscible in 1 mL of xylene  
Refractive index (n 20/D) : 1.515 - 1.525  
Density (at 25°C) : 0.915 - 0.925 g/mL  
Acidity : ≤ 0.05 ml

### Safety Information

Hazard Pictogram(s) :



Signal Word : Danger  
Hazard Statement(s) : EUH380- EUH430- H225- H304- H315- H319- H335-  
H336- H360FD- H373- H412  
Precautionary Statement(s) : P210- P260- P280- P301 + P310- P331- P370 + P378  
UN No. : 1993  
Class : 3  
Packing Group : II  
RTECS : Not available  
WGK : 3