

---

## Calcium bromide hydrate, Hi-LR™

**GRM6755**

### Product Identifier

---

CAS No.	:	71626-99-8
EC No.	:	232-164-6
Molecular Formula	:	CaBr <sub>2</sub> .xH <sub>2</sub> O
Molecular Weight	:	199.89 (anhydrous basis)
HS Code	:	2827 59 90
Storage	:	Below 30°C
Shelf life	:	4 years


### Technical Specification

---

Appearance	:	White crystals or powder or chunks or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (as CaBr <sub>2</sub> , EDTA Titration)	:	min. 83.00%

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P264- P280- P305+P351+P338- P337 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3