

Decanoic acid, For synthesis

GRM6758


Product Identifier

CAS No.	:	334-48-5
EC No.	:	206-376-4
Molecular Formula	:	C ₁₀ H ₂₀ O ₂
Molecular Weight	:	172.26
Synonym	:	Capric acid
HS Code	:	2915 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to pale yellow low-melting solid or liquid
Solubility	:	100 mg soluble in 1 mL of methanol
Assay (NaOH Titration/GC)	:	98.50 - 104.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H412
Precautionary Statement(s)	:	P264.- P273- P280- P302 + P352- P305+P351+P338- P332 + P313
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	HD9100000
WGK	:	1