
Zinc oxide, Hi-LRTM

GRM679


Product Identifier

CAS No.	:	1314-13-2
EC No.	:	215-222-5
Molecular Formula	:	ZnO
Molecular Weight	:	81.39
HS Code	:	2817 00 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White amorphous powder or solid
Solubility	:	33.3 mg soluble in 1 mL of dilute acids
Iron (Fe)	:	<= 0.04%
Assay (EDTA Titration)	:	98.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H410
Precautionary Statement(s)	:	P273- P391- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	ZH4810000
WGK	:	2