

Chromium sulphate basic

GRM6825


Product Identifier

CAS No.	:	39380-78-4
EC No.	:	235-595-8
Molecular Formula	:	$\text{Cr}_4(\text{SO}_4)_5(\text{OH})_2$
Molecular Weight	:	722.31
HS Code	:	2833 29 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Green to dark green powder
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Cr, ICP)	:	20.50 - 21.50%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P501
UN No.	:	1759
Class	:	8
Packing Group	:	II
RTECS	:	GB6240000
WGK	:	2