

---

## Molybdic acid, Hi-LR<sup>TM</sup>

**GRM689**

---

### Product Identifier

CAS No.	:	7782-91-4
EC No.	:	231-970-5
HS Code	:	2825 70 20
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to yellow crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL ammonium hydroxide solution
Chloride (Cl)	:	<= 0.01%
Sulfate (SO <sub>4</sub> )	:	<= 0.2%
Assay (MoO <sub>3</sub> , Lead nitrate Titration)	:	min. 85.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1

---

### Other Information

Additional information : Contains ammonium molybdate