
Molybdic acid, Hi-ARTM/ACS

GRM690

Product Identifier

CAS No.	:	7782-91-4
EC No.	:	231-970-5
HS Code	:	2825 70 20
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to yellow crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL ammonium hydroxide solution
Chloride (Cl)	:	<= 0.002%
Arsenate, phosphate and silicate (as SiO ₂)	:	<= 0.001%
Phosphate (PO ₄)	:	<= 0.0005%
Insoluble in dilute NH ₄ OH solution	:	<= 0.01%
Sulfate (SO ₄)	:	<= 0.2%
Heavy metals (as Pb)	:	<= 0.003%
Assay (MoO ₃ , Lead nitrate Titration)	:	min. 85.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1

Other Information

Additional information : Contains ammonium molybdate