

## 1,4-Dichlorobenzene

**GRM6907**


### Product Identifier

CAS No.	:	106-46-7
EC No.	:	203-400-5
Molecular Formula	:	C <sub>6</sub> H <sub>4</sub> Cl <sub>2</sub>
Molecular Weight	:	147.00
HS Code	:	2903 91 30
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White to off-white crystals or powder or solid or flakes
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Melting range	:	52 - 57°C
Assay (GC)	:	99.00 - 101.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H351- H410
Precautionary Statement(s)	:	P201- P273- P281- P308 + P313- P337 + P313- P501
UN No.	:	3077
Class	:	9
Packing Group	:	III
RTECS	:	CZ4550000
WGK	:	2