

GRM700 Potassium metabisulphite**Product Number Packing**

GRM700 : 500G

Product Information

Product Code : GRM700
Product Name : Potassium metabisulphite
Synonym : Potassium disulphite
Molecular Formula : $K_2S_2O_5$
Molecular Weight : 222.32
CAS No. : 16731-55-8
EC No. : 240-795-3
HS Code : 2832 20 00
Shelf Life : 4 years

Technical Specification

Appearance : White to light yellow crystals or powder
Solubility : 33.3 mg soluble in 1 mL of water
Assay (Iodometry) : 95.00 - 102.00%

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 31-36/37/38
S-Phrase(s) : 26-36
WGK : 1
RTECS : TT4920000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods