

Fast blue B salt

GRM7003


Product Identifier

CAS No.	:	14263-94-6
EC No.	:	238-153-2
C.I. No.	:	37235
Molecular Formula	:	C ₁₄ H ₁₂ N ₄ O ₂ Cl ₂ .ZnCl ₂
Molecular Weight	:	475.46
HS Code	:	3212 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to green to brown to black powder
Solubility	:	10 ppm in water yields clear solution
Absorption maxima 1	:	302 - 312 nm (10 ppm in water)
Absorption maxima 2	:	366 - 371 nm (10 ppm in water)

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H242- H350
Precautionary Statement(s)	:	P201- P210- P280- P308 + P313- P370 + P378- P403
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3