

Sodium metabisulphite, Hi-LR[®]

GRM707


Product Identifier

CAS No.	:	7681-57-4
EC No.	:	231-673-0
EC Index No	:	016-063-00-2
Molecular Formula	:	Na ₂ S ₂ O ₅
Molecular Weight	:	190.11
Synonym	:	Sodium pyrosulphite, Hi-LR [®]
HS Code	:	2832 10 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder or solid or chunks
Solubility	:	100 mg soluble in 1 mL of water
Iron (Fe)	:	<= 0.002%
Assay (Iodometry)	:	96.00 - 103.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302- H318
Precautionary Statement(s)	:	P270- P280- P305+P351+P338- P310- P330- P501
EU Phrases	:	EUH031
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UX8225000
WGK	:	1