

---

## Hydroxyethyl cellulose

**GRM7126**

---

### Product Identifier

CAS No.	:	9004-62-0
HS Code	:	3912 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White to light brown powder or powder with lumps
Solubility	:	33.3 mg soluble in 1 mL of water to give clear to hazy solution
pH (1% in water at 25°C)	:	6.00 - 8.50
Loss on drying (at 105°C, 3 hr)	:	<= 10.00%
Viscosity (5% in water at 25°C)	:	75 - 150 cps (Spindle No.: S1, 30 RPM)

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	1