
Hyflo[®] Super Cel[®]

GRM7136

Product Identifier

CAS No. : 68855-54-9
EC No. : 272-489-0
HS Code : 3802 90 19
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : White to light gray or pale buff powder
Solubility : Insoluble in water
pH (10% slurry in water at 25°C) : 8.5 - 11.0
Loss on ignition (at 750°C) : ≤ 2.00%

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1