

## Lauric acid, Hi-LR<sup>™</sup> For synthesis

**GRM7187**


### Product Identifier

CAS No.	:	143-07-7
EC No.	:	205-582-1
Molecular Formula	:	C <sub>12</sub> H <sub>24</sub> O <sub>2</sub>
Molecular Weight	:	200.32
HS Code	:	2915 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	White crystals or powder or flakes
Solubility	:	33.3 mg soluble in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Assay (NaOH Titration/GC)	:	99.00 - 102.00%

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H318
Precautionary Statement(s)	:	P280- P305+P351+P338+P310
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OE9800000
WGK	:	1