

Copper (II) carbonate basic

GRM721


Product Identifier

CAS No.	:	12069-69-1
EC No.	:	235-113-6
Molecular Formula	:	CuCO ₃ .Cu(OH) ₂
Molecular Weight	:	221.12
HS Code	:	2836 99 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Light blue to green crystal or powder
Solubility	:	100 mg soluble in 1 mL of 5N hydrochloric acid
Assay (Iodometry)	:	95.00 - 102.00%

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GL6910000
WGK	:	3