

---

## Magnesium stearate precipitated fine powder

**GRM7219**

### Product Identifier

---

CAS No.	:	557-04-0
EC No.	:	209-150-3
Molecular Formula	:	C <sub>36</sub> H <sub>70</sub> MgO <sub>4</sub>
Molecular Weight	:	591.24
HS Code	:	2915 70 90
Storage	:	Below 30°C
Shelf life	:	4 years


### Technical Specification

---

Appearance	:	White to off-white powder
Solubility	:	Insoluble in water
Assay (Mg, EDTA Titration)	:	3.80 - 5.00%

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	WI4390000
WGK	:	NWG