

---

## Manganese dioxide, Hi-AR<sup>®</sup>

**GRM7221**

### Product Identifier

---

CAS No.	:	1313-13-9
EC No.	:	215-202-6
Molecular Formula	:	MnO <sub>2</sub>
Molecular Weight	:	86.94
HS Code	:	2820 10 00
Storage	:	Below 30°C
Shelf life	:	4 years


### Technical Specification

---

Appearance	:	Grey to black crystals or powder
Solubility	:	Insoluble in water
Assay (KMnO <sub>4</sub> Titration)	:	99.00 - 102.00%

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302+H332
Precautionary Statement(s)	:	P261- P264- P270- P301 + P312- P304 + P340- P312
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OP0350000
WGK	:	2