

## Lead oxide, red

**GRM725**

### Product Identifier

CAS No.	:	1314-41-6
EC No.	:	215-235-6
Molecular Formula	:	Pb <sub>3</sub> O <sub>4</sub>
Molecular Weight	:	685.60
Synonym	:	Red lead oxide
HS Code	:	2824 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Red to orange powder
Solubility	:	Insoluble in water
Assay (EDTA Titration)	:	min. 95.00%

### Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H272- H302+H332- H351- H360- H372- H410
Precautionary Statement(s)	:	P201- P210- P260- P273- P280- P370 + P378
UN No.	:	1479
Class	:	5.1
Packing Group	:	II
RTECS	:	OG5425000
WGK	:	3