
Methyl cellulose (Low viscosity)

GRM7261

Product Identifier

CAS No. : 9004-67-5
HS Code : 3912 39 22
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : White to yellowish-white to greyish-white granules or powder or solid
Solubility : 33.3 mg soluble in 1 mL of water
Loss on drying (at 110°C, 1 hr) : ≤ 10.0%
Viscosity : 350 - 550 cps (c = 2% in water at 20°C)

Safety Information

UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : No data available
WGK : 1