

---

## Plaster of paris, Practical grade

**GRM7394**

---

### Product Identifier

CAS No.	:	10034-76-1
EC No.	:	231-900-3
Molecular Formula	:	CaSO <sub>4</sub> .½H <sub>2</sub> O
Molecular Weight	:	145.15
HS Code	:	2520 20 90
Storage	:	Below 30°C
Shelf life	:	4 years

---

### Technical Specification

Appearance	:	White with blue, grey or reddish tinge or brick red hygroscopic crystals or powder
Solubility	:	Insoluble in water
Full set time for water : powder (6:10)	:	~30 minutes

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1