

---

## Isopropyl myristate, Hi-LR™

**GRM7423**

### Product Identifier

---

CAS No.	:	110-27-0
EC No.	:	203-751-4
Molecular Formula	:	C <sub>17</sub> H <sub>34</sub> O <sub>2</sub>
Molecular Weight	:	270.45
HS Code	:	2915 90 99
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless to pale yellow liquid
Solubility	:	1 mL miscible in 1 mL of alcohol
FTIR	:	Matches with the standard pattern
Refractive index (n <sub>20/D</sub> )	:	1.4290 - 1.4390
Density (at 25°C)	:	0.845 - 0.855 g/mL
Assay (GC)	:	min. 98.00%

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	XB8600000
WGK	:	NWG