

GRM7466 Saccharin sodium dihydrate**Product Number Packing**

GRM7466 : 500G

Product Information

Product Code : GRM7466
Product Name : Saccharin sodium dihydrate
Molecular Formula : $C_7H_4NNaO_3S \cdot 2H_2O$
Molecular Weight : 241.19
CAS No. : 6155-57-3
EC No. : 204-886-1
HS Code : 2925 11 00
Shelf Life : 4 years

Technical Specification

Appearance : Colourless or white crystals or powder
Solubility : 100 mg soluble in 1 mL of water
FTIR (KBr disc) : Matches with the standard pattern
Loss on drying (at 120°C, 4hr) : $\leq 15.0\%$
Assay (NT) : 98.00 - 101.00%

Risk and Safety Information

WGK : nwg
RTECS : DE4375000
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : -
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods