

Salicylamide, Hi-LR™

GRM7469


Product Identifier

| | | |
|-------------------|---|---|
| CAS No. | : | 65-45-2 |
| EC No. | : | 200-609-3 |
| Molecular Formula | : | C ₇ H ₇ NO ₂ |
| Molecular Weight | : | 137.14 |
| HS Code | : | 2924 29 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------|---|---|
| Appearance | : | White to pink crystals or powder or solid |
| Solubility | : | 33.3 mg soluble in 1 mL of ethanol |
| FTIR | : | Matches with the standard pattern |
| Melting range | : | 139 - 143°C |
| Assay (NaOH Titration/GC) | : | min. 99.00% |

Safety Information

| | | |
|----------------------------|---|--|
| Hazard Pictogram(s) | : |  |
| Signal Word | : | Warning |
| Hazard Statement(s) | : | H302- H315- H319- H335 |
| Precautionary Statement(s) | : | P261- P305+P351+P338 |
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | VN6475000 |
| WGK | : | 3 |