
Shellac flakes, From Lacca

GRM7471

Product Identifier

CAS No.	:	9000-59-3
EC No.	:	232-549-9
HS Code	:	1301 20 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow to dark brown flakes
Solubility	:	Insoluble in water
Melting range	:	115 - 120°C
Saponification value	:	185 - 210

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3